



# Optec 5000P & 5500P

## Portable Vision Tester

- Ideal for industrial and DSE screening needs
- Dual lens system for testing bifocal patients
- 6/3 (20/15) to 6/60 (20/200) visual acuity range
- Unique maintenance free LED lighting system
- Portable
- No calibration required

The **Optec 5000P** and **5500P** are ideal for all industrial and DSE vision screening requirements. Visual acuity, depth perception, colour, phoria and peripheral vision can all be determined in a very efficient manner.

### Ease of use

A large viewing port, combined with a test pointer, facilitates fast and effective patient testing. Height adjustment ensures patient comfort and a remote control panel (**Optec 5500P**) saves time and gives excellent test flexibility.

### Accuracy

Several features provide for an extremely accurate vision test system. An advanced lighting system ensures high resolution, true colour test images with no glare reflection and there are no light bulbs to change as the instruments utilize a microprocessor controlled LED lighting system. A dual lens system permits a more accurate test on bifocal/multifocal patients.

### Comprehensive

Standardised test slide kits are available to meet the needs of industrial and primary care organisations. For specialist needs, there is a library of over 150 test slides available.

